



Clinical trial study summary: The DYNAMIC study


## Circulating tumor DNA guiding adjuvant therapy in stage II colon cancer<sup>1,2</sup>

Clinicians report the 2- and 5-year results from the DYNAMIC<sup>a</sup> study demonstrating the clinical benefit of MRD testing to inform therapeutic decisions for patients likely to benefit from adjuvant therapy.

### Background


In patients with stage II colon cancer, the benefit of adjuvant therapy is unclear, and identifying which patients should receive it is challenging. Some clinical studies have reported promising observational results on circulating tumor DNA (ctDNA) as a biomarker to measure minimal residual disease (MRD) postsurgery. The goal of the DYNAMIC study was to determine whether ctDNA-based MRD testing could improve the identification of 1) patients to receive adjuvant therapy and 2) patients who might forgo adjuvant therapy with minimal risk of recurrence. The DYNAMIC trial is the first to demonstrate clear clinical benefit of ctDNA-based MRD detection following surgery.

#### DYNAMIC study design



\*SOC = Standard of care

#### Key findings



### Study conclusions

Results from the study demonstrated that applying a ctDNA-guided approach to stage II colon cancer treatment management significantly reduced the use of adjuvant chemotherapy without compromising recurrence or survival at 2 years. Five-year follow-up results confirmed these findings, which also showed excellent survival outcomes with this approach.

<sup>a</sup> Haystack MRD™ uses an optimized version of the ctDNA detection technology used in the DYNAMIC study, which was the first study of its kind to assess the clinical benefit of MRD testing to guide adjuvant therapy. DYNAMIC is not a Haystack MRD-sponsored study.



<sup>1</sup> Tie J, Cohen JD, Lahouel K, et al. Circulating tumor DNA analysis guiding adjuvant therapy in stage II colon cancer. *N Engl J Med*. 2022;386(24):2261-2272. doi:10.1056/NEJMoa2200075

<sup>2</sup> Tie J, Wang Y, Serigne N L, et al. Circulating tumor DNA analysis guiding adjuvant therapy in stage II colon cancer: Overall survival and updated 5-year results from the randomized DYNAMIC trial. *J Clin Oncol*. 2024;42(6):supplement. doi:10.1200/JCO.2024.42.16\_suppl.108

# Gain *insights* that can help guide treatment options with more certainty than ever before

## Haystack MRD™

Haystack MRD™ is purpose-built for ultrasensitive ctDNA detection in patients with solid tumors and can be used to detect residual disease to guide adjuvant therapy decisions, to assess treatment response, and to monitor for recurrence.



## Better sensitivity means fewer false negatives and earlier insights

Haystack MRD enables ultrasensitive ctDNA testing with the ability to detect 95% of cases at 0.0005% tumor fraction, and in many cases, as low as 0.00008% tumor fraction, or approximately 1 ctDNA molecule in >1 million normal DNA molecules. With an exquisitely low limit of detection, insights can be uncovered to aid in earlier and more confident treatment decisions.

## The power of Quest Diagnostics

From early detection screening to diagnosis, monitoring, and beyond, Quest's oncology tests provide critical insights that help power personalized care at every step, making the transition between Haystack MRD and Quest's other advanced oncology testing seamless for easier, more streamlined patient care.



Contact our dedicated support team for Haystack MRD by emailing [providers@haystackmrd.com](mailto:providers@haystackmrd.com) or calling **1.844.966.7050**.

This test was developed and its performance characteristics determined by the CLIA-certified Haystack Oncology™ laboratory. It has not been cleared or approved by the US Food and Drug Administration.

Quest®, Quest Diagnostics®, any associated logos, and all associated Quest Diagnostics registered or unregistered trademarks are the property of Quest Diagnostics. © 2024 Quest Diagnostics Incorporated. All rights reserved. SB12657HAY 9/2024

